



## Laboratory diagnosis of hemostasis disorders

Diagnostic kits for the detection of SNP  
polymorphisms in genes associated with  
hemocoagulation and folate cycle



## Single Nucleotide Polymorphisms (SNP)



**Genetic polymorphism** – presence of at least two sequence variants (alleles) in a population, with a minor allele frequency of at least 1 %.

**SNP** – polymorphic traits represented by a single nucleotide substitution at a specific point in the genome. They are not an unconditional cause of the development of the disease, but can cause a greater or lesser risk of its development under the influence of various factors: SNP of other genes, concomitant diseases, pregnancy, lifestyle, smoking. It has been shown that SNPs contribute to the development of multifactorial diseases such as coronary heart disease, arterial hypertension and others.

### Multifactorial diseases

- determined by genetic and environmental factors
- most often they are polygenic, i.e. conditioned by SNP of several genes
- the detection of genetic risk factors indicates only a predisposition to a certain pathology, and not its presence
- environmental factors can be modified, and changing them can reduce the risk of developing the disease

One of the most studied genetic risk factors for developing multifactorial diseases are SNP genes associated with hemocoagulation and folate cycle.



### Indications for molecular-genetic testing for the presence of SNPs

- family history: family members with thrombotic and/or thromboembolytic complications, especially under 50 years of age,
- complicated obstetric history: difficulties during pregnancy, thrombotic difficulties after pregnancy and childbirth, two or more interruptions in fetal development in the early stages, birth of children with malformations, etc.,
- planned preparation for pregnancy or extracorporeal fertilization (especially in women who have a history of thrombosis or a family history of thromboembolic complications)
- planned hormonal contraceptives or hormone replacement therapy,
- planned massive surgical interventions, long-term immobilization.

## RealBest® Technology: Solution for differential diagnosis of SNP in genes of hemocoagulation system and folate cycle.



### Biological specimen

- **whole blood**
- **buccal epithelium** (buccal epithelium collection is collected into the Transport Solution (Cat. No. 8894))



### Preparation

- Automatic: **RealBest GenMag** (Cat. No. 8846)  
Manual: **RealBest-Genetics DNA-express** (Cat. No. 8845)



### Amplification

- 🧪 **Detection of 2 SNP in a single tube:**
  - RealBest-Genetics Hemostasis (F2/F5)**
  - RealBest-Genetics Hemostasis (FGB/F13A1)**
  - RealBest-Genetics Hemostasis PAI-1/ITGB3**
  - RealBest-Genetics Hemostasis ITGA2/F7**
- 🧪 **Detection of 4 SNP in two tubes:**
  - RealBest-Genetics Hemostasis (MTR/MTRR/MTHFR)**
- 🧪 **Detection of 12 SNP in four tubes:**
  - RealBest-Genetics Hemostasis (12)**



### Interpretation

Genotype identification: **healthy homozygote/heterozygote/mutant** homozygote is performed automatically.

## Basic characteristics of SNP in genes of hemocoagulation system and folate cycle

Gene/SNP	Mutation effect	Aggravating factors
<b>MTHFR</b> (methylenetetrahydrofolate reductase), <b>C677T</b>	10x higher risk of developing hyperhomocysteinemia for homozygotes (T/T). Risk associated with the development of thrombosis and atherosclerosis, pregnancy difficulties and fetal malformations	<i>F5</i> : G1691A; <i>F2</i> : G20210A; Deficiency of vitamins B6, B12, folic acid, defects of other genes of the folate cycle; smoking
<b>MTHFR A1298C</b>	Decrease in enzyme activity, the risk of developing hyperhomocysteinemia	
<b>MTR</b> (methionin synthase), <b>A2756G</b>	Decrease in methionine concentration and accumulation of homocysteine in the bloodstream: increased risk of cardiovascular disease, pregnancy problems	<i>MTRR</i> : A66G; <i>F5</i> : G1691A; <i>F2</i> : G20210A; deficiency of B6, B12, folic acid; smoking
<b>MTRR</b> (methionine synthase reductase), <b>A66G</b>	Decrease in enzyme activity, increased risk of developing hyperhomocysteinemia. The risk of developing thrombotic complaints, fetal malformations	<i>MTR</i> : A2756G; <i>F5</i> : G1691A; <i>F2</i> : G20210A; deficiency of B6, B12, folic acid; smoking
<b>F2</b> ( <i>coagulation factor II</i> , <i>prothrombin</i> ), <b>G20210A</b>	Increased level of prothrombin in plasma: increased risk of developing venous thrombosis, pregnancy difficulties	<i>F5</i> : G1691A; smoking
<b>F5</b> ( <i>coagulation factor V</i> , <i>proaccelerin</i> ), <b>G1691A - Leiden mutation</b>	Increased risk of thrombosis and embolism: heterozygotes (G/A) have a 7x higher risk, homozygotes (A/A) have a 20x higher risk. Risk of developing birth and gynecological problems	Hormonal contraceptives, protein C and S deficiency hyperhomocysteinemia; <i>F2</i> : G20210A; smoking
<b>F7</b> ( <i>coagulation factor VII</i> , <i>proconvertin</i> ), <b>G10976A</b>	Reduction of factor VII concentration in plasma, 2x reduction in the risk of developing thrombosis and myocardial infarction even in coronary atherosclerosis	Absence of <i>F5</i> : G1691A; <i>F2</i> : G20210A; Presence of <i>F13</i> : G103T
<b>F13</b> ( <i>fibrin-stabilizing factor XIII</i> ), <b>G103T</b>	Formation of thinner fibrin clots: reduced risk of developing cardiovascular diseases (thrombosis, myocardial infarction)	Absence of <i>F5</i> : G1691A; <i>F2</i> : G20210A; <i>FGB</i> : G(-455)A; Presence of <i>F7</i> : G10976A
<b>PAI-1</b> ( <i>plasminogen activator inhibitor-1</i> ), <b>5G(-675)4G</b>	Decrease in fibrinolytic activity of blood: increased risk of thrombosis, as well as the development of obstetric and gynecological problems	<i>ITGB3</i> : T1565C; <i>F5</i> : G1691A; <i>F2</i> : G20210A; ACE: del1287
<b>ITGA2</b> ( <i>integrin alpha-2</i> ), <b>C807T (c.C759T)</b>	Increased platelet adhesion: homozygotes (T/T) have been shown to be at increased risk of developing thrombotic complications	<i>F5</i> : G1691A; <i>F2</i> : G20210A
<b>ITGB3</b> ( <i>integrin beta-3</i> ), <b>T1565C (c.T176C)</b>	Platelet hyperaggregation: increased risk of developing thrombosis and embolism. Increased risk of rupture of atherosclerotic plaques. In carriers of the T allele, the effectiveness of using aspirin is reduced.	<i>F5</i> : G1691A; <i>F2</i> : G20210A
<b>FGB</b> ( <i>coagulation factor I</i> , <i>fibrinogen</i> ), <b>G(-455)A</b>	Increase in plasma fibrinogen concentration: increased risk of developing ischemic and cardiovascular diseases, difficulties during pregnancy and the postpartum period	<i>FGB</i> : G(-455)A; smoking

## Diagnostic kits for the detection of SNP polymorphisms in genes associated with hemocoagulation and folate cycle

Cat. №	Kit name	Number of tests	Gene: Polymorphisms
<b>Extraction kits for the isolation of nucleic acids</b>			
8845	RealBest-Genetics DNA-express	50	-
8846	RealBest GenMag	96 (2x48)	-
<b>SNP detection kits</b>			
3802 €€	RealBest-Genetics Hemostasis (F2/F5)	48	F2: 20210G/A; F5: 1691G/A
3803 €€	RealBest-Genetics Hemostasis (MTR/MTRR/MTHFR)	48	MTR: 2756A/G; MTRR: 66A/G; MTHFR: 677C/T; MTHFR: 1298A/C
3831 €€	RealBest-Genetics Hemostasis (FGB/F13A1)	48	FGB: (-455)G/A; F13A1: c.103G/T
3832 €€	RealBest-Genetics Hemostasis ITGA2/F7	48	ITGA2: 807C/T; F7: 10976G/A
3833 €€	RealBest-Genetics Hemostasis PAI-1/ITGB3	48	PAI-1: -6755G/4G; ITGB3: 1565T/C
3801	RealBest-Genetics Hemostasis (12)	48	F2: 20210G/A; F5: 1691G/A; F7: 10976G/A; F13A1: c.103G/T; FGB: (-455)G/A; ITGA2: 807C/T; ITGB3: 1565T/C; PAI-1: -6755G/4G; MTR: 2756A/G; MTRR: 66A/G; MTHFR: 677C/T; MTHFR: 1298A/C

### Features and advantages of the diagnostic kit

- **Ready Master Mix for PCR:** Simplification of analysis procedures and high stability of test quality
- **Universal protocol:** Multiple different tests in one run
- **Multiplexes:** determination of two SNPs in one tube
- **Specimens:** whole blood or buccal epithelium
- **Compatible devices:** CFX96 (Bio-Rad, USA), Gentier 96E/R (Xi'an TianLong, Science and Technology Co., Ltd., China)
- **High stability of the kit:** storage at a temperature of 2–8 ° C; transport up to 26 ° C not more than 10 days

Diagnostické centrum DNK, s.r.o.

Brestová 14, 821 02 Bratislava

+421 911 299 324, +421 911 211 404

[dnk@pharma.sk](mailto:dnk@pharma.sk), [diagnostika@pharma.sk](mailto:diagnostika@pharma.sk)

[www.pcr.sk](http://www.pcr.sk)